




Economical, Reusable & Moisture Resistant Plastic

True North® Film freezer boxes by Heathrow Scientific, were designed with sustainability in mind and are made from a specially formulated thin film polypropylene. True North freezer boxes are moisture resistant, which significantly reduces the chance of mold growth. With superior durability, thin film polypropylene will outlast cardboard boxes that are susceptible to softening and producing dust and provide a vastly more cost effective alternative to polycarbonate boxes. Safe to use down to -196°C (vapour phase), they are compatible with low temperature storage vessels and mechanical freezers. True North Freezer Boxes are compatible with most standard 2" (50mm) freezer racking.

When it comes to saving valuable lab space, True North Freezer Boxes are supplied "flat-packed", which eliminates the hassle of stacking and reduces the amount of space needed for storage. When ready to use, assembly is easy and pre-assembled dividers are included with each box. True North freezer boxes come complete in packs of 10, and are available in four different colors.

Thin Polypropylene film 1.5-2 mL

Part No.	Tube	Wells	Well Dia.	L x W x H	Pkg.	Color	
HS120370	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	White	□
HS120371	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	Blue	■
HS120372	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	Lt. Green	■
HS120373	1.5/2.0 mL	81	12.1 mm	4.95 x 4.92 x 1.92 in (12.5 x 12.4 x 4.9 cm)	10/pk	Purple	■



Freezer boxes are supplied "flat-packed"

Can be used in mechanical freezers and vessels down to -196°C (Vapour Phase only)

