



# SPZV-50 Stereo Zoom Trincoular Microscope

with Mighty Cam Pro on Double Arm Boom Stand

This high definition 1080P auto focus inspection system with a 22" monitor offers superb image quality at 60 frames per second. You can operate the system with or without a PC.

This innovative all-in-one inspection system gives operators two viewing options:

#### SPZV-50 Trinocular Microscope

- View objects beneath the microscope in crystal-clear 3D, with a magnification range of 6.7x - 50x

### Mighty Cam Pro - HDMI Monitor or PC

- The Mighty Cam Pro inspection camera connects directly to your HD monitor via HDMI, or to your PC via USB
- View real-time images at a magnification range of 34x 255x

### **Key Features**

- Stereo Zoom Technology Provides a crystal-clear 3D view of objects below the microscope
- Dual View- Examine objects on the HD monitor or through the microscope eye-pieces
- HDMI Mode Save images, measure, and more without a PC
- SD Card Included Save image and video files in HDMI mode
- PC Mode Connect to a PC to utilize more advanced software features

Dual View Inspect through the monitor or eye-pieces HDMI & USB Connect to an HD monitor or PC SD Card Included

Save images/ video files Double Arm Boom Stand Easily move across the object being



## HDMI Mode - Integrated Measurement Software No PC Required



## USB Mode - Capture Software Connect to a PC for more advanced features



**Vertical and Horizontal Adjustment** 

Stereo Zoom Trinocular Microscope SPZV-50 [6.7x - 50x] with Mighty Cam Pro on Double Arm Boom Stand	
Item #	26800B-383-PRO
Magnification	Microscope Magnification Range: 6.7x - 50x On-Screen Magnification Range: 34x - 255x
Resolution	Dynamic Resolution: 2MP 1920 x 1080P Static Resolution: 6MP 3264 x 1836P
Illumination	Integrated LED Ring Light
Output Interface	<ul><li>USB</li><li>HDMI</li><li>SD Card</li></ul>
Package Contents	Mighty Cam Pro (26100-258)  Stereo Zoom Microscope (SPZV-50)  Double Arm Boom Stand (26800B-534)  E-Arm LED Focus Mount with brightness control box (26200B-209)  22" LCD HD Monitor (26700-406)